

Search for Transglutaminase 2 (TG2) covalent inhibitors using Domainex's covalent fragment library, and a LC-MS based covalent fragment screening platform

Scott Martin, Anna Hopkins, Gary Wu, Jason Pembroke, Megan Ashley and Andrew Ratcliffe

Domainex Ltd, Chesterford Research Park, Little Chesterford, Saffron Walden, CB10 1XL

Introduction


Recent advances in the tools to support discovery/design of covalent inhibitor drugs and the success of block buster drugs such as Osimertinib and Ibrutinib has led to increased interest in 'electrophile-first' covalent drug discovery. Covalent fragment-based screening by intact protein mass spectrometry (MS) has been shown to be a powerful tool, demonstrated by the discovery of KRAS(G12C) inhibitors [1]. Additional assays to support covalent fragment-based screening, including assessing warhead reactivity via a GSH assay and binding site identification by proteolytic digestion and peptide mapping further optimise hits. Covalent inhibition is time-dependent, so the preferred measure of potency is the second-order rate constant k_{inact}/K_i , rather than IC50.

Here we describe a case study assessing covalent fragments against Transglutaminase 2 (TG2) protein (a protein that is implicated in the pathogenesis of several diseases, including cancer, fibrosis and neurodegenerative diseases), which to our knowledge has not been done before.

Covalent Screening by MS


Initial screening of pooled fragments

- Following method optimisation, TG2 was incubated with pools of fragments for 18 hours at +4°C and covalent binding was measured by rapid intact LC-MS under denatured condition.
- Fragment pools were prepared using an automatic sample handler, fragments were pooled into groups of 5 with >5 Da mass difference.
- LC-MS data were acquired on Waters G2-XS QToF, utilising the chromatography from a Waters ACQUITY UPLC Protein BEH C4 VanGuard Pre-column, 300Å, 1.7 µm, 2.1 mm X 5 mm on a Waters Acquity H-Class Plus Bio. Data of up to 800 pools can be acquired per night.
- LC-MS data were analysed and visualized by an in-house developed data processing pipeline, which can process 4 000 LC- HRAMS data files within 20 minutes.



Selecting individual fragments to validate hits

- Example data of fragment bindings in individual runs vs pooled runs
- Side-by-side data comparisons were performed by an in-house developed rapid data processing pipeline.




Glutathione Counter screen

- Hits from the covalent screen were assessed for Glutathione (GSH) reactivity in real time by LC-UV-MS
- Compounds with high reactivity ($T_{1/2} < 100$ minutes) are unlikely to exhibit desirable properties and were discarded at this stage.

Compound concentration	1.67 µM
Final solvent concentration	<0.2% DMSO
Positive control	N-nitrophenyl acrylamide
Incubation duration	8 hours*
Analysis	LCMS with UV-DAD
Compound requirements	50 µL of 10 mM DMSO stock or 1 mg of solid material
Deliverables	half-life ($T_{1/2}$), presence of GSH adduct (Y/N), PBS stability (%)

Fragment	$T_{1/2}$	PBS Stability (% remaining)
A	>5000	110
B	>5000	117
C	2923	109
D	>5000	110
E	>5000	114
F	1023	114
G	1118	113
H	148	110
I	>5000	106
J	473	114

- Multiple fragments show relatively high binding with low GSH reactivity
- All 10 strong binder hits had a $1/2$ life >100 minutes and were taken forward for binding site identification by proteolytic digestion, followed by peptide mapping.
- Exemplar data for fragment A is shown below




Binding Site Identification

- Fragments were individually incubated with TG2 for 18 hours at +4°C
- The samples were then denatured with urea, reduced with TCEP, alkylated with IAA and digested overnight at +37°C with trypsin

Digestion of the protein/fragment complex into peptides

LC-HRAMS analysis




Peptide mapping software identifies the peptide where fragment A is bound

```

AHHHHHAAE  LVLRCDLEL  ETNGRDHHTA  DLCREKLVVR  RGQPFWLTLH  FEGRNYEASV  DSLTFSVVTG  PAPSQEAGTK
EGDWTATVVD  QDCTLSLQL  TTPANAPIGL  YRLSLEASTG  YQGSFVLGH  FILLFNANCP  ADAVYLDSEE  ERQEYVLTQQ
IKNIPHNFQ  FEDGILDICL  ILLDVNPKFL  KNAGRDCSRR  SSPVYGRVV  SGMVNCDDQ  GVLGRLWNN  YDGVSPMSH
KNHGCRVKY  GQCVFAAVA  CTVLRCLGIP  TRVVTNNSA  HDQNSLLIE  YFRNFEQIQ  GDKSEMIWF  HCWVESMTR
QALDPTPQEK  SEGTYCCGPV  PVRAIKEGDL  STKYDAPFV  AEVNADVVD  IQQDGSVHK  SINRSLVGL  KISTKSVGRD
YPEGSEERE  AFTRANHLK  L
    
```


MS/MS data is used by the software to identify the exact binding site within the peptide



- Binding site data was successfully generated for all 10 fragments
- All fragments were shown to bind to the target cysteine (residue 277) of TG2

K_{inact}/K_i

- Covalent inhibition is time-dependent, so the preferred measure of potency is the second-order rate constant k_{inact}/K_i , rather than IC50
- TG2 was incubated with a fragment at different concentrations for different periods of time. The binding percentage at each time point was measured by LC-MS.
- An in-house fast data processing pipeline has been developed to automatically plot the binding data and extract a k_{inact}/K_i value for each compound



Conclusions

- Domainex has developed a MS based covalent platform which successfully identified several fragments as potential leads for optimisation against TG2.
- GSH counter screen of positive hits can be used to triage out highly reactive warheads
- Proteolytic digestion and peptide mapping determine the binding site of selected hits
- k_{inact}/K_i was successfully determined for selected hits and can be used in combination with GSH reactivity to optimise compound potency during the hit to lead stage
- Applications of these workflows allows for the high throughput screening of various libraries against a target protein, as well as supporting decisions during of the optimisation process.

Future Development

- Evaluate a LC-MS based covalent fragment screening platform for alternative binding sites, beyond cysteine.

Contact

If you would like to learn more about our drug discovery platforms, please contact: enquiries@domainex.co.uk